## **DICOM CD/DVD startup guide (ENGLISH)**

This short guide explains how to launch the viewer application which allows displaying the DICOM medical images present on this CD/DVD. Depending on the configuration of the application which created this CD/DVD, the viewer application may also be absent from the present media.

Instructions for each supported operating system follow.

### ***Windows operating system***

If your PC runs a Windows operating system, the *LocalEye Viewer* application may start automatically once the CD/DVD is inserted in the drive. If this doesn't happen (for instance, due to the fact that autorunning of removable media is disabled by the Operating System for security reasons), please browse to the root of CD/DVD contents, and double-click on the "*ViewCDWindows.bat*" file: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

### ***Mac OS X operating system***

If you are using a Mac, please insert the CD/DVD in the appropriate drive, wait a few seconds, then double click on the CD/DVD icon which will appear on your desktop, in order to open the

CD/DVD and browse its contents. Double-click on the "*ViewCDMacOSX.command*" file located in the root of the CD/DVD: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

### ***Linux operating system***

If you are using a Linux PC with a graphic windowing environment, please insert the CD/DVD in the appropriate drive, wait a few seconds, then double click on the CD/DVD icon which will appear on your desktop (auto-mounting shall be enabled), in order to open the CD/DVD and browse its contents. Double-click on the "*ViewCDLinux.sh*" file located in the root of the CD/DVD: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

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